

	No	Type		No	Type		No	Type
AB Dental			Bio3			Cowellmedi		
Narrow Platform	55021	11	Advanced	55067	59	Inno Int Octa	55014	4
Standard Platform	55049	41	Progressive	55067	59	Inno Sub Hex	55037	29
ADIN			BioHorizons			Dentack		
Touareg S/OS	55034	26	Laser-Lok 3.0 Grey	55034	26	Internal Hex	55056	48
Touareg Closefit UNP 2.75	55059	51	Internal Laser-Lok 3.5 Yellow	55049	41			
Touareg Closefit NP 3.0	55059	51	Internal Laser-Lok 4.5 Green	55051	43	Dental Ratio		
Touareg Closefit RP 3.5	55046	38	Internal Laser-Lok 5.7 Blue	55051	43	Tissue Level RP, WP	55014	4
Touareg CloseFit WP 4.3, 5.0	55067	59	External 4.0 Platform	55011	1	Bone Level NC	55061	53
Swell 3.3	55034	26				Bone Level RC	55062	54
AIRES			Biomet/3i			Dental Tech		
3.3	55049	41	External Hex	55011	1	Implassic FT3 3.25	55039	31
			Certain internal Purple	55036	28	Implassic FT3 3.75	55049	41
			Certain internal Purple Platform 3.4	55036	28			
			Certain internal Blue, Yellow, Green	55066	58	Dentalis		
			T3 Purple	55036	28	Regular	55051	43
			T3 Blue, Yellow, Green	55066	58			
Alfa Gate			Bionnovation			Dentegris		
Bioactive, All diam	55034	26	Biomorse	55023	14	3.3	55034	26
MAX, All diam	55034	26				3.75	55057	49
M+ NP 3.5	55044	36				4.5	55049	41
M+ RP 4.3 - 5.0	55037	29				3.8/4.5	55049	41
Porous, All diam	55034	26				5.5	55020	10
Alpha Dent			Biotech Dental			Dentis		
Conical connection Narrow Platform	55025	16	Contact, Contact-S 3.6 - 4.8	55045	37	Cleanlant S-Clean 3.7 - 4.8	55020	10
Conical connection Wide Platform	55057	49				S-Clean OneQ Narrow 3.3	55064	56
Hexagon connection Narrow Platform	55039	31	BlueSkyBio					
Hexagon connection Regular Platform	55057	49	BIO Max NP Platform, All	55046	38	Dentium		
			BIO Max RP Platform 4.3-5.0	55037	29	Superline	55037	29
			BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	NR Line	55066	58
Alpha-Bio			BIO Int Hex 3.5 Platform 4.7	55051	43	Simple Line II GH 1,6	55014	4
Internal Hex	55034	26	BIO Int Hex 4.5 Platform 4.7-7.0	55014	4			
Conical Hex 3.2	55065	57	BIO One-stage Narrow, Regular, Wide	55014	4	DEVA		
Conical Hex 3.75-5.0	55051	43	BIO One-stage Short Regular & Wide	55014	4	Bone Level 3.5	55064	56
			BIO Quattro Narrow 3.3	55061	53	Bone Level 4.5	55051	43
Anker			BIO Quattro Regular 4.1 - 5.6	55062	54			
SB III 3.5	55023	14	BIO Conus 12, 3.0 Platform 3.0	55064	56	DIO		
SB III 4.0-5.0	55020	10	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38	SM Submerged 3.8	55049	41
			BIO Conus 12, 4.5 x 6mm	55046	38	SM Submerged 5.0	55020	10
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7	UF II Narrow	55064	56
			BIO Conus 12, 5.0 x 6mm	55046	38	UF II Regular	55020	10
			BIO Trilobe 3.5	55016	6			
Ankylos			BIO Trilobe 4.3-6.0	55017	7	Ditron Dental		
C/X	55025	16				MPI 3.3-6.0	55087	79
			Bredent			Ultimate 3.75-6.0	55087	79
			SKY	55025	16			
Anthogyr			copaSKY	55027	19	Easy Implant		
Axiom Bone Level 3.4 - 5.2	55065	57	mini2SKY	55033	25	Ocean	55046	38
Axiom Tri-lobe	55065	57	fast & fixed Abutment	55083	75	Lilas	55037	29
Axiom Tissue Level Narrow, Regular	55050	42						
			BTI			Edierre		
Argon			Interna Universal	55051	43	3-lobe	55021	11
K3pro Rapid 3.0 - 3.5	55026	17	Interna Wide	55051	43	4-lobe	55025	16
K3pro Rapid Platf 3mm	55019	9	BTK					
K3pro Rapid Platf 3mm Short	55044	36	External Hex EN	55038	30	Euroteknika		
			External Hex ER	55011	1	Internal Platform 3.0	55064	56
			External Hex EW	55011	1	Internal Platform 3.5	55081	73
			Internal Hex IR	55049	41	Naturactis/Natural+ 3.0	55064	56
			Internal Hex IM	55051	43	Naturactis 3.5-5.0	55081	73
			Internal Hex IW	55051	43	Uneva	55011	1
			Hexagonal Morse Taper KR	55081	73			
Astra Tech			Hexagonal Morse Taper KW	55056	48	galimplant		
Osseospeed EV 3.0	55058	50	Octagonal Morse Taper SR	55014	4	iPx 4012	55046	38
Osseospeed EV 3.6	55046	38						
Osseospeed EV 4.2	55057	49	C-Tech			GC Tech		
Osseospeed EV 4.8	55055	47	EL	55067	59	Aadva Narrow	55081	73
Osseospeed EV 5.4	55048	40				Aadva Regular	55066	58
Osseospeed TX 3.0 S	55064	56	Camlog					
Osseospeed TX 3.5-4.0 S	55046	38	Screw-Line, Root-Line 3.3	55036	28	Geass		
Osseospeed TX 4.5-5.0 S	55017	7	Screw-Line, Root-Line 3.8, 4.3	55031	23	Way Extra	55081	73
			Screw-Line, Root-Line 5.0, 6.0	55032	24			
			iSy	55065	57	Global D		
Avinent						In-Kone 4.0-5.0	55024	15
External, Coral & Ocean 3.5	55038	30	Ceralog					
External, Coral & Ocean 4.1	55011	1	Hexalobe ceramic	55076	68			
External, Coral 5.1	55013	3						
Internal, Coral & Ocean 3.5, 4.1	55051	43	Conelog					
Conical, Ocean 3.5	55066	58	Screw-Line 3.3	55066	58			
Conical, Ocean 4.1	55004	18	Screw-Line 3.8, 4.3	55066	58			
			Screw-Line 5.0	55075	67			
B&B Implant								
Duravit Slim	55021	11						
Duravit 3P	55035	27						
Duravit Evolution	55035	27						
BEGO								
Semados Mini	55012	2						
Semados S / RI 3.25-5.5	55034	26						
Semados RS / RSX 3.0	55082	74						
Semados RS / RSX 3.75-5.5	55034	26						
Semados SC / SCX 3.25-5.5	55034	26						

	No	Type		No	Type		No	Type
gmi			Implant Swiss			Medentika		
avantgard RP 3.75-4.25	55066	58	Tissue level	55014	4	Microcone NI	55021	11
avantgard WP 4.75-5.75	55075	67	Bone level Regular	55037	29	Microcone RI	55081	73
frontier RP,WP	55051	43	Bone level 3.3	55065	57	Quattrocone RI	55081	73
phoenix NP 3.3	55038	30				Quattrocone 30 AI	55081	73
phoenix Std 3.3 - Wide 5.0	55011	1	Implantvision			Procone 3.3	55119	85
Monolith 3.0	55043	35	Internal Hex 3.7-5.7	55034	26	Procone 3.8	55031	23
						Procone 4.3	55031	23
						Procone 5.0	55032	24
GP Implants			Innobiosurg			Medentis		
SSI Spiral Shape implant	55051	43	MagiCore Regular	55018	8	ICX Premium 3.75 - 4.8	55062	54
						ICX Plus 3.45	55039	31
						ICX Narrow	55058	50
gt medical			Intoss Anchor			Medical Instinct		
HXE	55011	1	Internal Hex 3.7-5.7	55088	80	Hex 3.0	55071	63
HXE RP	55012	2				Hex 3.4, 4.0, 5.0	55052	44
HXE WP	55013	3	Intra-lock			Cone 3.0, 3.4, 4.0	55071	63
HXI 4.1	55037	29	Intrahex 3.4 - 4.5	55051	43	Cone 5.0	55072	64
002 HXI 4.7	55049	41	Gold & Blue 3.4 - 6.0	55051	43			
OCT 4.3	55014	4	InDex 3.4 - 3.75	55015	5			
			InDex 4.0 - 4.3	55020	10			
			FlatTop 3.75	55012	2			
Hahn			IRES			Medigma Biomedical		
Tapered 3.0	55058	50	External Hex Narrow 3.3	55012	2	Mars 3,5-6,0	55057	49
Tapered 3.5	55046	38	External Hex Regular 3.7 - 5.2	55011	1	FixTite 3,5-6	55057	49
Tapered 4.3	55046	38	Internal Hex Narrow 3.3	55087	79	Fix-a-Dent 3,75-6	55057	49
Tapered 5.0	55037	29	Internal Hex Regular 3.7 - 5.2	55049	41			
Tapered 7.0	55048	40	Internal Octagon Regular 3.7 - 5.2	55014	4	Medimecca		
						PT SLA Chaorum	55047	39
Hi-Tec			J Dental Care			Megagen		
Tapered Self Thread NTI 3.3, SLT 3,75	55034	26	JD Evolution 3.7-6.0	55067	59	Anyridge	55051	43
Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	JD Evolution Plus+ 3.7-6.0	55067	59	Anyone Internal	55037	29
Spark TPI	55034	26	JD Evolution S 3.25	55044	36	Mini	55064	56
Logic Plus	55034	26	JD ICON ultra, S 2.75	55064	56	MUA	55080	72
Logic LGI 3.0-3.5	55046	38	Conical Abutment	55080	72			
Logic LGI 4.3-5.0	55037	29						
Hiossen			Keystone			Microdent		
HS System	55028	20	Prima / Connex	55049	41	Universal	55011	1
ET System	55037	29	Prima / Connex	55067	59	Microdent 2.8 - 3.25	55068	60
ET 3.2	55047	39	Prima / Connex	55041	33	Microdent 3.3 - 5.0	55011	1
ET mini	55015	5	Prima / Genesis Narrow	55049	41	Microdent 5.5 - 8.0	55013	3
			Prima / Genesis Regular	55067	59	Genius 3.0	55025	16
IDI Evolution			Prima / Genesis Wide	55041	33	Genius 3.5 - 5.0	55051	43
Standard implant 3.75-4.0	55004	18	Restore Narrow	55079	71	Ektos 3.2	55023	14
B1One High	55080	72	Restore Regular	55011	1	Ektos 3.7 - 8.0	55051	43
B1One Low	55033	25	Restore Wide	55013	3	Trylogic 3.5 - 4.2	55051	43
						Trylogic 5.0	55075	67
IDI Implants			Klockner			MIS		
ID Cam ST 4.2	55010	13	VEGA Micro 3.0	55074	66	Conical connection Narrow (C1, V3)	55081	73
ID All	55010	13	VEGA Narrow 3.5	55065	57	Conical connection Standard (C1, V3)	55087	79
ID3 3.1	55021	11	VEGA Regular 4.0, 4.5	55034	26	Conical connection Wide (C1)	55041	33
			KL Narrow 3.3	55012	2	Seven and M4 Narrow	55027	19
ids			KL Regular 3.7, 4.2	55011	1	Seven and M4 Standard	55049	41
Reflect Aspire	55081	73	KL Wide 4.7	55013	3	Seven Wide 5.0	55051	43
Reflect Certus	55019	9	SK2	55038	30	Connect abutment NP, SP, WP	55115	81
Reflect Tapered	55034	26	Essential Cone	55067	59			
Reflect Recover	55020	10	S3M 3.1	55079	71			
Reflect Rapid	55020	10				Neo Biotech		
			Konmet			IS II 4.0	55020	10
Implance			Classic 4.0, 4.8	55014	4	IS III 4.0	55020	10
Tissue level	55014	4	Cylindrical 3.5-4.0	55065	57			
Bone level Regular	55037	29	Cylindrical 4.5-5.5	55057	49	Neodent		
Bone level 3.3	55065	57	Short Cylindrical 5.0-6.0	55078	70	Drive CM	55025	16
			Conical 4.4/3.4	55065	57	Titamax WS	55056	48
Implant Direct			Conical 4.8/3.8, 5.5/4.5	55057	49	GM	55066	58
Legacy 1/2/3/4 3.0/3.2	55012	2	Kyocera			Neoss		
Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43	Bone Level 3.2-3.4 (NP)	55058	50	SP (Ø 3.5 implants and wider)	55004	18
RePlant 3.5	55022	12	Bone Level 3.7-4.2 (RP)	55034	26	NP (Ø 3.25)	55027	19
RePlant 4.3 - 5.0	55010	13	Tissue Level 3.7-4.7 (RP)	55014	4	Access Abutment	55060	52
RePlant 6.0	55032	24	Tissue Level 4.7-5.2 (WP)	55014	4			
ReActive 3.5 Platform	55022	12				Nobel Biocare		
ReActive 4.3 - 5.0 Platform	55010	13	Lasak			Brånemark System Ext hex RP	55011	1
Interactive 3.0 Platform	55046	38	BioniQ QR	55066	58	Brånemark System Ext hex NP	55012	2
Interactive 3.4 Platform	55037	29	BioniQ QN	55021	11	Brånemark System Ext hex WP	55013	3
Screwindirect 5.0 PF 3.0 - 5.7	55070	62				Tri-lobe Connection NP	55022	12
			Leader			Tri-lobe Connection RP, WP	55010	13
Implant One			Internal Hex, 3.3 - 5.5 mm	55034	26	Tri-lobe Connection 6.0	55032	24
100 Series Std thread	55029	21	Cono-In	55067	59	Conical Connection 3.0	55058	50
200 Series Std thread	55029	21				Conical Connection Narrow	55046	38
300 Series Std thread (blue)	55029	21	Maxillent			Conical Connection Regular	55037	29
400 Series Std thread (purple)	55035	27	iRaise	55049	41	Conical Connection Wide 5.5	55048	40
500 series Std thread (green)	55035	27				MUA NP & RP	55080	72
300 Series Wide thread (blue)	55029	21				On1 NP, RP, WP	55042	34
400 Series Wide thread (purple)	55035	27						
500 Series Wide thread (green)	55035	27						

	No	Type		No	Type		No	Type
Noris Medical			P-I Brånemark			Spiral Tech		
Tuff 3.3 - 6.0	55034	26	HEX 3.3	55115	81	ESi, Hex connection	55051	43
			HEX 3.75 -5.0	55117	83	Ultimate, Hex connection	55051	43
Novem			AMP 3.3	55039	31	ESi Conical 3.5	55044	36
3.4	55039	31	AMP 3.75 - 4.8	55118	84	ESi Conical 4.3, 5.0	55037	29
3.8	55016	6	MT 3.3 - 4.8	55046	38			
4.2, 5.0	55057	49	Conical Abutment HEX, AMP, MT	55080	72	Straumann		
						Tissue level, Regular & Wide Neck	55014	4
NucleOss			Prodent Italia			Tissue level, Narrow Neck Crossfit (NNC)	55076	68
Tpure Bone Level	55004	18	Twinner & Prime, Conomet, Collar: Orange	55027	19	Bone level / Bone Level Tapered Narrow	55061	53
T6 Bonel Level	55074	66	Twinner & Prime, Conomet, Collar:: Fuchsia	55051	43	Bone level / Bone Level Tapered Regular	55062	54
			Twinner & Prime, Conomet, Collar: Green	55067	59	Bone level / Bone Level Tapered SC 2.9	55064	56
OCO Biomedical			Twinner & Prime, Conomet, Collar: Yellow	55020	10	Bone Level X (BLX) all platform	55027	19
Engage 3.25	55027	19	Prime, Conomet, Collar: Blue	55020	10	Screw retained Abutment NP, RP	55080	72
Engage 4.0-6.0	55051	43						
ERI 3.25	55069	61	Quest			Surgikor		
ERI 4.0-5.0	55014	4	Hylock NP	55030	22	Fixation 3.0	55058	50
TSI 3.25	55068	60	Hylock RP	55066	58	Fixation 3.5	55066	58
TSI 4.0-5.0	55011	1				Fixation 4.3	55037	29
			Resista			Solution 3.5, 4.0	55066	58
OsseoFuse			WIS 3.0-5.0	55049	41	Solution 4.5, 6.0	55037	29
Hexaplus/ HexaPlus-S Sinus 3.75-7.5	55081	73				Immediate, all	55049	41
ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5	Rex Implants			Versatile, all	55049	41
			External hex TL 1.8	55060	52			
Osstem			Internal hex TL 2.9	55034	26	Sweden & Martina		
US Regular	55011	1				Syra, Syra SP 3.9-5.0	55011	1
US Mini	55012	2	Ritter			Outlink ² 3.75, 4.10, 5.0	55011	1
US Wide	55013	3	Multi Unit	55033	25	Shelta 3.8-5.0	55025	16
US Wide PS	55011	1	Spiral SB/LA Narrow PF 2.9	55081	73	Premium Straight & SP 3.3-3.8	55025	16
SS Regular, Wide	55028	20	Spiral SB/LA Standard PF 3.75	55049	41	Kohno Straight & SP 3.3-3.8	55025	16
TS Regular, Reg Short, Diam 6 & 7	55037	29				Premium Straight & SP 4.25-6.0	55004	18
TS Mini 3.5	55015	5	RS Implant			Kohno Straight & SP 4.25-6.0	55004	18
TS Mini 3.0	55047	39	3.5 - 6.0	55020	10	Prama 3.8 - 5.0	55079	71
						DAT Connection	55035	27
Osteo Plus			Schutz Dental			Tag Dental		
Iridium PF 3.8 & 4.2	55079	71	Impla cylindrical	55052	44	Axis	55051	43
She 3.45, 3.75	55011	1	IMPLA micro Retention Cone 4.2	55052	44	massif	55051	43
Shi 3.75, 4.2	55049	41						
Shi 5.0	55051	43	SGS Implant			Tatum		
			P1, P7 Platform	55049	41	Integrity 3.7 - 5.0	55015	5
OT-Medical			P1D, P7D	55049	41	Integrity 6.0 - 8.0	55040	32
OT-F2 3.4, 3.8	55053	45	P7N	55059	51	T implant	55063	55
OT-F2 4.1	55073	65						
			SIC invent			TBR		
Overmed			MUA	55120	86	Infinity PF 3.5, PF 4.0	55004	18
Easyline Easydip	55020	10	SiCace 3.4 - 5.0	55052	44	M Implant	55044	36
ETLINE	55034	26	SiCmax, SiCtapered 3.7 - 5.2	55052	44			
			SiCvantage® max, SiCvantage tapered 3.0	55026	17	Thommen		
Oxy Implant			SiCvantage® max, SiCvantage tapered 3.7	55116	82	SPI 3.5	55021	11
Internal	55020	10	SiCvantage® max, SiCvantage tapered 4.2-5.2	55054	46	SPI 4.0	55065	57
Kone Conical Narrow	55021	11				SPI 4.5-6.0	55066	58
Kone Conical Regular	55051	43	Simply Implants					
PSK Line Regular	55051	43	Internal Hex	55051	43	ticare		
External Mini	55079	71				Osseus Mini, Standard	55011	1
External Regular	55011	1	SIN			Osseus Maxi	55013	3
External Maxi	55011	1	Strong SW CM	55023	14	Osseus Maxi ps	55011	1
			Strong SW Int Hex	55035	27	Inhex Mini	55021	11
Paltop			Unitite CM 2.9	55084	76	Inhex Standard	55046	38
Standard 3.75-5.0	55034	26	Unitite CM 3.5	55077	69	Inhex Maxi	55017	7
Narrow	55027	19	Unitite CM Compact 4.0	55085	77	Osseus Quattro Standard	55011	1
PCA Conical	55066	58	Tryon 3.25 - 5.0	55011	1	Inhex Quattro Standard	55046	38
Phibo			Southern Implants			TRATE		
Aurea EVO NP	55027	19	Tri-Nex 4.3-5.0	55010	13	ROOTFORM	55065	57
Aurea EVO RP, WP	55051	43	Tri-Nex 3.5	55022	12			
TSA S3	55074	66	External Hex	55011	1	TRI Dental implants		
TSA S4, S5	55034	26	MAX Ext Hex	55011	1	Tissue level Octa	55014	4
TSH S3, S4, S5	55011	1	Co-Axis Ext Hex	55011	1	Tri-Vent, all	55034	26
			Piccolo Ext Hex	55050	42	Narrow	55025	16
			IT-Connection	55014	4			
			Deep Conical 3.0	55064	56			
			Deep Conical 3.5-4.0	55046	38			
			Deep Conical Tapered 5.0	55017	7			

	No	Type
Trinon		
Q2-Implant	55011	1
QK-Implant	55014	4
QZA 3.5 Platform	55049	41
QZA 4.5, 5.7 Platform	55051	43
URIS		
Omni Narrow 3.0-3.5	55115	81
Omni Regular 4.0-6.5	55020	10
Vulkan		
Internal Hex	55034	26
External Hex NP	55012	2
External Hex RP	55011	1
Conical con MiP	55058	50
Conical con NP	55046	38
Winsix		
Internal Hex	55020	10
Xgate		
Conical 3.3 - 5.0	55034	26
Xive		
S-Plus 3.0	55030	22
S-Plus 3.4	55061	53
S-Plus 3.8-5.5	55062	54
TG Plus 3.4-4.5	55018	8
Zeramex		
XT Implant 4.2, 5.5	55086	78
Ziacom		
Conical connection RP	55046	38
Conical connection WP	55075	67
External hexagonal connection NP	55079	71
External hexagonal connection RP	55117	83
External hexagonal connection WP	55014	4
Internal hexagonal connection NP	55027	19
Internal hexagonal connection RP	55057	49
Internal hexagonal connection WP	55057	49
Zimmer		
Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43
Tapered Swiss Platf 3.8	55070	62
Tapered Swiss Platf 4.8	55014	4
Spline Twist Platf 4 & 5	55038	30
Eztetic NP 3.1	55065	57
Zuga		
All	55027	19